

NEXT-GENERATION SELF-ADHESIVE SUPERABSORBER

curea P1 easy border

silicone-free wound contact surface · multipurpose




cureamedical

Innovation in Care

curea P1 easy border – adheres securely, supports cleverly

curea P1 easy border is a self-adhesive multifunctional wound dressing. The silicone adhesive border securely and gently fixes the wound dressing, while the silicone-free wound contact surface can absorb exudate of various viscosities without obstruction.



Completely reliable and comfortable – for worry-free wound care and a better quality of life.



Size	Max. wound size	Item No. 10 pcs
12 x 12 cm	5.5 x 5.5 cm	P1EB-120120-10

More sizes on the way.



Secure hold

curea P1 easy border stays securely in place, even with movement. Reliable adhesion you can trust.



Gentle on skin

The silicone adhesive border can be removed easily and atraumatically during dressing changes, while the integrated notches (FLEXI CUTS) ensure a comfortable wearing experience.



Reliable exudate management without compromise

curea P1 easy border quickly absorbs and safely retains all types of exudate (even under compression), helping to reduce the risk of maceration.



Easy handling

curea P1 easy border is self-explanatory to apply and can be repositioned if necessary. No additional fixation material is required.




Versatile in use

curea P1 easy border is suitable for light to heavy exuding wounds, offering a wide range of applications.

Convincing in quality and performance – layer by layer

User- and patient-friendly design with FLEXI CUTS allows for easy and flexible application, and is also comfortable to wear and less constricting.

FLEXI CUTS



PU-backing film

Breathable • Liquid-impermeable • Easy application • Elastic and flexible

Breathable Textile Backsheet (BTBS)

Sterile barrier against germs from inside and outside • Protection against leakage and contamination

Soft absorbent core

Low to high exudate levels • Exudate of varying viscosities • Reliable absorption and secure retention • Well-balanced moist wound healing environment • Free of adhesives and binders

Silicone-free nonwoven wound contact layer

Highly permeable • No silicone on the wound • Unhindered exudate absorption without barriers • Supports wound bed cleansing (soft debridement)

Surrounding silicone adhesive border

Secure hold, but gentle on the skin • No additional fixing material required • Easy to remove and reposition

Typical areas of application

- ✓ **Pressure ulcer**
- ✓ **Leg ulcer**
- ✓ **Postoperative wounds**
- ✓ **Skin Tears**



Next-generation superabsorber

Like all our wound dressings, curea P1 easy border also bears the 'MPD' seal. MPD dressings represent the next generation of superabsorbents, supporting simple, comprehensive and economical care of exuding wounds as all-in-one products.

